

# JUST ADD T.L.C.

To give your respiratory patient the very best protection against auto-infection

The Control III® Disinfectant Home Care Kit is so complete, all you have to add is T.L.C.

Our double-quat disinfectant is completely effective against germs in the hospital and in home respiratory therapy. And it's odorless, economical and easy to use!

All your patient has to do is dilute the Control III Disinfectant concentrate with water and soak pre-cleaned equipment for 10 minutes. The solution can be stored and reused for up to 14 days. Test strips are provided to check solution activity.

Just add your special brand of T.L.C. to give your patients everything they need for successful home care.



## **CONTROL III® Disinfectant**

CONTROL III Disinfectant is a time-proven germicide concentrate formulated to disinfect equipment in the home. CONTROL III is a double quaternary ammonium chloride compound, containing a special formulation of germ killing compounds effective against a broad spectrum of bacteria, viruses and fungi.

### **Available in a Home Care Kit**

- 2 ounce bottle of CONTROL III Disinfectant (concentrate)
- 2-quart container with lid
- Measuring cup
- Test Strips & color key
- Instructions

### **Time Proven**

- Used since 1974

### **Easy to Use**

- Works in 10 minutes
- Odorless
- Non-irritating
- Non-staining

### **Economical**

- Up to 8 times less expensive than vinegar
- Lasts up to 14 days
- 2-quart container can be reused

### **Effective**

- Highly effective against gram-negative and gram positive bacteria, viruses and fungi
- Effective against Pseudomonas, Vaccina, HIV-1, Influenza A<sub>2</sub>, Herpes Simplex 1, Escherichia, Staphylococcus, and more
- Residual bacteriostatic film
- EPA approved as hospital disinfectant

**CONTROL III** concentrate is also available in 8 and 16 ounce bottles, and in ready-to-use gallons. 3 and 5.5 gallon soak buckets and test strips are available separately.

Put T.L.C. into your patient care. Contact Maril Products and order your Home Care Kits today.



## **MARIL PRODUCTS, INC.**

320 West 6th Street  
Tustin, CA 92780  
800-546-7711  
714-544-7711  
714-544-4830 (fax)

[control3@earthlink.net](mailto:control3@earthlink.net) (e-mail)